

## Nuaire LabGard 605 Small Animal Portable Isolator

The LabGard 605 Small Animal Portable Isolator provides protection for animals, technicians, or both simultaneously. The 605 creates a mini clean room to isolate infected or sterile product.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Nuaire](#)

### Description



Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

- [Features](#)
- [Specifications](#)
- [Optional Features](#)

## Features

### Recycled Airflow

The LabGard 605 creates a mini clean room through recirculated airflow. The recycling feature produces ultra clean air for both personnel and animals.

### HePeX™ Zero Leak Airflow System

The HePeX™ plenum inflates to provide quiet uniform air without transferring vibration. Positive pressure chambers and ducts are surrounded by negative pressure eliminating potential leaks. HePeX™ fabrics are incredibly durable. They are flameproof, provide a barrier against water, gases, fluid, and dust as well as resist to chemicals, oils, and greases. It even holds antibacterial properties and is impregnable meaning harmful materials become trapped in the filters.

Motor Speed Control

Supply HEPA filter 99.99% @ 0.3 microns

Prefilter

Protective Filter Grill

Stainless Steel Animal Rack

Stainless Steel Castors

4 Stainless Steel Animal Shelves

Charcoal Filter

Specifications

Model	NU-605-500	NU-605-600
Nominal Size	5 ft. (1.5 m)	6 ft. (1.8 m)
Style of Cabinet	Console	Console
Exterior Dimensions (W x D x H)	66 x 79 x 78	78 x 79 x 78
in. [mm]	(1676 x 2007 x 1981)	(1981 x 2007 x 1981)
Weight	750 lbs. / 340 kg.	800 lbs. / 398 kg.
Inflow Velocity	30 fpm (0.15 m/s)	30 fpm (0.15 m/s)
Warranty	3 Years Parts and Labor (US &CA) 4 Year Parts (Intl)	

Accolades UL, UL-C, CE, Personnel Protection, Product Protection, RoHS  
Please note specifications are subject to change. Please contact NuAire to verify before purchase.

Optional Features

- Optional Features**  
Additional Stainless Steel Rack  
Additional Stainless Steel Shelves  
Clean Room Canopy

Reviews

There are yet no reviews for this product.