

Galileo RapidCast 0911 Horizontal Gel Box

The RapiCast 0911 gel electrophoresis system has an extended run-length to improve resolution and separation of macromolecules.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Galileo](#)

Description

The RapiCast 0911 gel electrophoresis system has an extended run-length to improve resolution and separation of macromolecules. This gel box also features a RapiCast gasket tray for fast, leak-proof, in-unit gel casting.

System Description

The 0911 mini gel system includes a 9cm wide x 11cm long RapiCast™ gasketed gel tray for convenient gel casting right in the buffer chamber.

System features

- Heavy duty lower buffer chamber
- Interlock safety lid with attached leads
- Gasketed (cast in place) 9cm wide x 11cm long UV-transparent gel tray with built-in fluorescent ruler and 2 comb slots.
- 2 combs of your choice:

80-0911 Minigel System, gasketed UVT gel tray, 9cm x 11cm & (1) 10 and (1) 14 tooth comb, 1.5 mm T

80-0911- Comb 1.0mm thick/5 tooth

C5-1.0

80-0911- Comb 1.5 mm thick/5 tooth

C5-1.5

80-0911- Comb 1.0mm thick/8 tooth

C8-1.0

80-0911- Comb 1.5 mm thick/8 tooth

C8-1.5

80-0911- 1X Microtiter format Comb 1.0 mm thick/9 tooth

CMT9-1.

0

80-0911- 1X Microtiter format Comb 1.5 mm thick/9 tooth

CMT9-1.

5

80-0911- Comb 1.0 mm thick/10 tooth
C10-1.0
80-0911- Comb 1.5 mm thick/10 tooth
C10-1.5
80-0911- Comb 1.0 mm thick/12 tooth
C12-1.0
80-0911- Comb 1.5 mm thick/12 tooth
C12-1.5

Reviews

There are yet no reviews for this product.