

## Benchmark BV101 Vornado™ Miniature Vortex Mixer

Powerful vortexing for tubes up to 50mL. Smooth, quiet and efficient operation. Despite its powerful motor, the Vornado has a footprint of less than 4x4 inches, allowing it to fit on even the most crowded bench. Unique head design prevents liquids from entering the housing, prolonging motor life.

Rating: Not Rated Yet  
[Ask a question about this product](#)

Manufacturer[Benchmark Scientific](#)

### Description

Powerful is the best word to describe this mini vortexer. With a 4mm orbit and fixed speed of 2800rpm, the Vornado instantly vortexes even the largest samples, including nearly full 50 ml tubes. Despite its powerful motor, the Vornado has a footprint of less than 4x4 inches, allowing it to fit on even the most crowded bench.

Insta-Touch operation means that the motor is automatically activated with slight pressure on the center cup head. Operation is smooth, quiet and efficient. When pressure is removed, operation immediately ceases.

The Vornado's unique head design prevents liquids from entering the housing, prolonging motor life. The unit is available with five different cup head colors.

- Powerful vortexing for tubes up to 50mL
- Insta-Touch™ pressure activated operation
- Mini size - less than 4 inches wide
- Smooth, quiet and efficient

## Benchmark BV101 Vornado™ Miniature Vortex Mixer Specifications

Speed:	2800 rpm
Orbit:	4mm
Motor Activation:	Insta-Touch™ pressure activation
Dimensions: (w x d x h)	3.7 x 3.9 x 2.6 in. 9.4 x 9.9 x 6.6 cm
Weight:	1 lbs. / 0.4 kg
Electrical:	Universal Voltage 100-240V, 50-60 Hz, 0.7A
Warranty:	2 Years

### Reviews

There are yet no reviews for this product.