

## Helmer iLR111 Laboratory Refrigerator

The iLR111 is a single door i.Series® Laboratory Refrigerator. It features the i.C<sup>3</sup> full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description ??

- [Description](#)
- [Features](#)
- [Specification](#)

### Description

The iLR111 is a single door i.Series® Laboratory Refrigerator. It features the i.C<sup>3</sup> full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

### Features

- Innovative i.C<sup>3</sup> User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i Plus warranty includes 7 years compressor, 2 years parts, and 1 year labor

## Specification

Interior Dimensions (w x h x d)

- 20" x 49.5" x 20"
- 508mm x 1258mm x 508mm

Exterior Dimensions (w x h x d)

- 24.25" x 70.5" x 24.5"
- 616mm x 1791mm x 623mm

Interior Volume (w x h x d)

- 11.5 cu ft / 326 liters

Shelves (w x d)

- 4 Adjustable Full-Size
- 18" x 18" (458mm x 458mm)

Electrical

- 115V 60Hz 7 FLA
- 230V 50Hz 3.5 FLA
- 230V 60Hz 3.5 FLA

Net Weight

- 322 lb / 147 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

## Reviews

There are yet no reviews for this product.