

Helmer HLR104-ADA Undercounter Laboratory / Pharmacy Refrigerator

The HLR104-ADA is a single door ADA compatible Undercounter Horizon Series Laboratory/Pharmacy Refrigerator. It fits under 34" ADA compliant counter heights. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The HLR104-ADA is a single door ADA compatible Undercounter Horizon Series Laboratory/Pharmacy Refrigerator. It fits under 34" ADA compliant counter heights. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Features

- Temperature controller with alarm/monitor system provides protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Rel.i warranty includes 5 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 19.5" x 17.75" x 22"
- 496mm x 451mm x 559mm

Exterior Dimensions (w x h x d)

- 24" x 31.5" x 26.75"
- 610mm x 801mm x 680mm

Interior Volume (w x h x d)

- 4 cu ft / 113 liters

Shelves (w x d)

- 2 Adjustable Full-Size
- 18" x 18" (458mm x 458mm)

Electrical

- 115V 60Hz 5 FLA

Net Weight

- 185 lb / 84 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

Reviews

There are yet no reviews for this product.