

## Helmer HLF120 Laboratory Freezer

The HLF120 is a single door Horizon Series Laboratory Freezer. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

### Description

- [Description](#)
- [Features](#)
- [Specification](#)

## Description

The HLF120 is a single door Horizon Series Laboratory Freezer. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

## Features

- Temperature controller with alarm/monitor system provides protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Rel.i warranty includes 3 years compressor, 2 years parts, and 1 year labor

## Specification

Interior Dimensions (w x h x d)

- 24.75" x 58.25" x 24.25"
- 629mm x 1480mm x 616mm

## Single Door Upright: Helmer HLF120 Laboratory Freezer

---

### Exterior Dimensions (w x h x d)

- 29.5" x 80" x 29.5"
- 750mm x 2032mm x 750mm

### Interior Volume (w x h x d)

- 20.2 cu ft / 572 liters

### Shelves (w x d)

- 4 Adjustable Full-Size
- 22" x 23" (559mm x 585mm)

### Electrical

- 208 / 230V 60Hz 4.3 FLA
- 230V 50Hz 3.8 FLA

### Net Weight

- 440 lb / 200 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

### Reviews

There are yet no reviews for this product.