

Helmer iPR120 Pharmacy Refrigerator

The iPR120 is a single door i.Series® Pharmacy Refrigerator. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The iPR120 is a single door i.Series® Pharmacy Refrigerator. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Features

- Innovative i.C³ User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i Plus warranty includes 7 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 24.75" x 58.25" x 24.25"

Single Door Upright: Helmer iPR120 Pharmacy Refrigerator

- 629mm x 1480mm x 616mm

Exterior Dimensions (w x h x d)

- 29.5" x 79.5" x 29.5"
- 750mm x 2020mm x 750mm

Interior Volume (w x h x d)

- 20.2 cu ft / 572 liters

Shelves / Baskets (w x d)

- 1 Adjustable Full-Size
- 22" x 23" (559mm x 585mm)
- 6 Roll Out
- 21" x 21.25" (534mm x 540mm)

Electrical

- 115V 60Hz 7.5 FLA
- 230V 50Hz 4.2 FLA
- 230V 60Hz 4.2 FLA

Net Weight

- 525 lb / 239 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

Reviews

There are yet no reviews for this product.