

Steelco AC 7000 NHP

For large animal cages and racks. The NHP series is a heavy duty, large capacity, hydro-spray cage and rack washer (1200mm/47.24" chamber width available).

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Steelco](#)

Description

- [Description](#)
- [Specification](#)

Description

The AC 7000 NHP series cage and rack washer is a heavy duty, large capacity, hydro-spray washer specifically designed to support the throughput of medium and large size facilities with large animal cages and racks.

The machine is standard equipped with three preheating water tanks and two dosage units resulting in great cycle settings flexibility.

AC 7000 NHP is available in single or double interlocked hinged door made of stainless steel AISI 316L with glass inspection window. Door orientation can be chosen both right or left hand.

Useful chamber dimensions: Width 1200mm (47.24"), Depth 2240mm (88.19"), Height 2350mm (92.52")

Specification

- Automatic manifold for trolleys equipped with spray arms or injection systems (optional for feeding bottles treatment)
- Flush self-cleaning filter for the recirculation of the washing water.
- Double chemical control dosing system.
- Data storage system and process control.
- Operator multi language touch screen control with auto-diagnostic checking that constantly monitors and displays current cycles status and alarms, allows to optimize the washing process and to personalize programs on board of the machines. Up to 65 washing and rinsing programs available: 5 service programmes, 10 memorized programmes and 50 free programmes.

External dimensions: Width 2350mm (92.52"), Depth 2488mm (97.95"), Height 2720mm (107.08")

Reviews

There are yet no reviews for this product.