

Nuair LabGard ES 819 Class I Enclosure

The LabGard ES 819 Class I Enclosure is designed to contain and remove contaminants from the work surface with the aid of an external exhaust system for applications involving liquids, powders, animal necropsy, or equipment.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Nuair](#)

Description



Ergonomic Design Contributes to a friendly working environment



Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

- [Features](#)
- [Specifications](#)

Features

Choice in Airflow

The 819 can be used as a dead air box, connected for a facility exhaust system, or be used with NuAire's HEPA Filter Blower Module to create a Class I Biological Safety Cabinet

Two Service Ports
Hinged Viewing Acrylic Window
Clear Acrylic Side Walls
White Acrylic Back Panel
Coved Black Epoxy Black Work Surface

Optional Features

Low Flow Alarm
Base Support Stands
Exhaust HEPA Filter Blower Module
Flexible Exhaust Duct
Seismic Brackets

Specifications

Model

Nominal Size

Sash Height

Exterior Dimensions (W x D x H)
in. [mm]

Interior Dimensions (W x D x H)
in. [mm]

Weight

Warranty	3 Year Parts and Labor (US & CA) 4 Year Parts (Intl)
----------	---

Accolades	UL, CE, Energy Saver, Personnel Protection, Environmental Protection
-----------	--

Please note specifications are subject to change. Please contact NuAire to verify before purchase.

Reviews

There are yet no reviews for this product.