



Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rockers and MORE!



INCUSHAKER™ WITH

SHAKING INCUBATORS

MAGiC CLAMP™
PATENTED

- Exceptional temperature uniformity
- Large flask capacity
- MAGiC Clamp accessories for instant-magnetic clamp exchange
- Available with CO₂ control or cooling
- **NEW!** Adjustable clamp for 125, 250 & 500ml flasks



Item: H1000-M

INCUSHAKER
Mini

Speed Range: 30 to 300rpm
Temp Range: Amb. +5° to 70°C
Platform Dims: 9.5 x 11.5 in.
Universal Platform: H1000-MR
Flask Capacity: 2 x 1L (H1000-MR-1000)
5 x 500ml (H1000-MR-500)
8 x 250ml (H1000-MR-250)
12 x 125ml (H1000-MR-125)

\$2895



Item: H1010

INCUSHAKER
10L

Speed Range: 30 to 300rpm
Temp Range: Amb. +5° to 60°C
Platform Dims: 17.8 x 17.8 in.
Universal Platform: H1010-MR
Flask Capacity: 5 x 2L (H1000-MR-2000)
9 x 1L (H1000-MR-1000)
16 x 500ml (H1000-MR-500)
25 x 250ml (H1000-MR-250)
36 x 125ml (H1000-MR-125)

\$5675



Item: H1012

INCUSHAKER
10LR (with cooling)

Speed Range: 30 to 300rpm
Temp Range: Amb. -15° to 60°C
Platform Dims: 17.8 x 17.8 in.
Universal Platform: H1010-MR
Flask Capacity: 5 x 2L (H1000-MR-2000)
9 x 1L (H1000-MR-1000)
16 x 500ml (H1000-MR-500)
25 x 250ml (H1000-MR-250)
36 x 125ml (H1000-MR-125)

\$7995



Item: H3501

INCUSHAKER
Mini CO₂

Speed Range: 30 to 300rpm
Temp Range: Amb. +5° to 60°C
CO₂ Range: 0 to 20%
Platform Dims: 9.5 x 11.5 in.
Universal Platform: H1000-MR
Flask Capacity: 5 x 500ml (H1000-MR-500)
8 x 250ml (H1000-MR-250)
12 x 125ml (H1000-MR-125)

\$8495



Limited Time Offer*
FREE!
myTEMP^{Mini}

\$459
Value

Item: H2200-H
(Ambient +5° to 60°C)



*With purchase of any Incu-Shaker

Special offers are valid through July 31st 2018. Promo Code: BP2.18 must be referenced at time of order.



ESPIRE
LABORATORY ASPIRATOR

- Compact Design
- Built-in Vacuum Pump
- Adjustable Vacuum Pressure
- Vacuum Gauge
- Safe and Precise Aspiration

\$825

V0020	Aspire™ Laboratory Aspirator	\$1995.00
Prefilled kits: 50 x 2.0ml tubes, prefilled with high impact bead media		
V1002-1SL	Optional single channel adapter, extra long (80mm) stainless steel (D)	\$146.00
V1002-8	Optional eight channel adapter, for use with 200ul pipette tips	\$146.00
V1002-58	Optional eight channel adapter, fine tip stainless steel (E)	\$146.00

Limited Time Offer*
FREE!

\$359
Value

Combi-Rotor™
- 12 x 1.5ml
- 4 x PCR Strips

myFUGE™ 12
Mini Centrifuge
C1012

*With Purchase of
V0020 or B3D5000
Use Promo Code BP2.18



BENCHWAVER™
Undulating 3-D Shaker

- Undulating 3D motion
- Adjustable speed and tilt
- Ideal for staining, washing and general mixing



Optional Stacking Platform

\$1179

Speed: 5 - 105rpm
Platform: 13.5 x 13 in.
Tilt Angle: 0 to ±10°

B3D5000	BenchWaver™ Undulating Shaker	\$1179.00
B3D5000-STK	Optional stacking platform	\$119.00
B3D5000-DIMP	Optional dimpled mat	\$79.00

Also available . . .

BenchMixer™ V2



- Reversible Top:
Use "Cup" or "Flat/Dimpled" surface
- Powerful motor for INSTANT vortexing
- Variable speed with touch or continuous operation



BV1003
BV1000

BenchMixer™ V2 Vortex Mixer (with Reversible Top)
BenchMixer™ Vortex Mixer (with standard cup head)

\$265.00
\$239.00



Iso-BLOCK™

- Two isolated chambers
- Digital temperature control
- Hinged lid
- Available with Quick Flip™ universal blocks

Serves as two baths in one!

BSH6000
BSWCMB

IsoBlock™ digital dry bath, without blocks
Block, Quick Flip, 24 x 1.5/2.0ml, 14 x 0.5ml,
6 x PCR strips/ 48 x 0.2ml

\$679.00
\$129.00

*Please visit website for additional block availability.

ORBi-SHAKER™



19mm
Orbit

- Digital speed control, 30 - 300rpm
- High capacity (9x500ml or 14x250ml)
- Exchange flask clamps in 1 second
with MAGic Clamp™ system

With large platforms and a wide range
of accessories, the Orbi-Shaker is useful
for a variety of applications in the
general research laboratory.

BT3000
BT3000-MR

Orbi-Shaker with non-slip rubber mat, 12x12"
MAGic-Clamp platform for flasks & tubes, 14x12"

\$1494.00
\$229.00

*Additional accessories are available.

myBATH™ SERIES

- Digital accuracy and convenience
- Quick-CAL™ user calibration function
- Seamless, stainless steel chamber



Rack for
29mm tubes

Rack for
17mm tubes

Racks included with
4L water bath, optional
with 8L and 12L models



B2000-2
B2000-4
B2000-8
B2000-12

myBath 2 Liter water bath
myBath 4 Liter water bath
myBath 8L water bath
myBath 12L water bath

\$545.00
\$795.00
\$869.00
\$1039.00

LABORATORY BALANCES



Precision Balance

- Readability 0.001 g
- Easy in-lab calibration feature
- RS232 port for data output to a
computer or printer



Series Tx

- High Resolution Touch Screen
- USB port for flash drive data acquisition
- RS32 port for computer / printer connection
- Quick-Cal™ Internal Calibration System

W3200-120
W3200-320
W3200-500
W3002A-120
W3002A-220

Precision Balance, 120g x 0.001g
Precision Balance 320g x 0.001g
Precision Balance 500g x 0.001g
Analytical Balance series Tx, 120 x 0.0001g
Analytical Balances series Tx, 220g x 0.0001g

\$925.00
\$925.00
\$1075.00
\$2995.00
\$2295.00

NEW! MC-24™ Touch

High Speed Microcentrifuge



- Full Color Touch Screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 in.
- Variable speed control, set in rpm or rcf



Max. Speed: 13,500rpm / 16,800xg
Max Capacity: 24 x 1.5/2.0ml
2 x PCR Strips (16 x 0.2ml)

C2417 MC-24™ Touch Microcentrifuge
with 24 place COMBI-Rotor

\$1995.00

Digital Hotplates & Stirrers



- White, ceramic top plate
- 200 to 1500rpm stir speed
- Ambient+5°C to 380°C temp. range
- Optional temperature probe

These hotplates and stirrers feature a large
LED display and digital control for unmatched
accuracy and performance.

H3760-S
H3760-H
H3760-HS

Magnetic stirrer, 200-1500rpm
Hotplate, Ambient+5°C to 380°C
Hotplate with magnetic stirrer

\$219.00
\$219.00
\$339.00