



# Monitor and record

New Brunswick™ BioCommand® SFI Track and Trend

# BioCommand® SFI Track and Trend

PC-Based Monitoring and Data Logging for New Shakers, Freezers & CO<sub>2</sub> Incubators



Now you can monitor and permanently record the key operating parameters from multiple shakers, CO<sub>2</sub> incubators, freezers, or other OPC-compliant systems throughout your lab from a single PC.

BioCommand® SFI Track and Trend software is an easy-to-use, inexpensive, PC-based software that can trend and archive critical data including temperature, shaking speed, gas concentrations and humidity for up to 32 units, simultaneously.

## BioCommand® SFI Advantages

- > Provides full monitoring and data logging capability
- > Enables rapid setup through intuitive programming
- > Offers text alarm notifications
- > Allows alarms to be annotated for complete traceability using a user-defined "add event" feature
- > Lets you customize trend screens, to monitor eight critical parameters of your choosing, per window
- > Compatible with Windows® 7, Windows Vista® and Windows XP® operating systems

Catalog No.	BioCommand SFI Packages*
M1291-1000	Shaker Package
M1291-1001	CO <sub>2</sub> Incubator Package
M1291-1002	ULT Freezer Package

\* Includes software, and cable to connect one shaker, freezer or CO<sub>2</sub> incubator. Computer requires an RS-232 or USB port. Shaker/incubator/freezer requires RS-232 or RS-485 port.

Catalog No.	To Connect Additional Units**
M1291-8001	Shaker Cable, 25 ft./7.6 m
P0620-7012	CO <sub>2</sub> Incubator Cable, 15 ft./4.5 m
M1291-5001 †	ULT Freezer Cable Set, 45 ft./15 m
P0460-7751	4-Port RS-232/485-to-USB Converter
P0460-7750	8-Port RS-232/485-to-USB Converter

\*\* When connecting multiple shakers, CO<sub>2</sub> incubators or freezers, in addition to a single software package, one additional cable is required for each additional system you monitor/log. Shakers and CO<sub>2</sub> incubators additionally require a 4- or 8-port converter. † Includes RS-485 Board.

Your local distributor: [www.eppendorf.com/worldwide](http://www.eppendorf.com/worldwide)  
 Eppendorf AG · Hamburg · Germany · Tel: +49 40 538 01-0  
 Application Support E-mail: [support@eppendorf.com](mailto:support@eppendorf.com)  
 Eppendorf North America, Inc. · USA · Tel: +1 516 334 7500  
 Toll free phone: +1 800 645 3050 menu option 2  
 E-mail: [techserv@eppendorf.com](mailto:techserv@eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany. New Brunswick™ is a trademark of Eppendorf AG, Germany. BioCommand® is a registered trademark of New Brunswick Scientific Co., Inc. Windows®, Windows Vista® and Windows XP® are registered trademarks of Microsoft Corporation in the United States and other countries. All rights reserved, including graphics and photos. Order no B-4030. Copyright © 2013 by Eppendorf AG.