

forte

A.T. VILLA

adaptable laboratory system



Attractive
Affordable
Adaptable

A.T. Villa's Forte system stands uncontested as one of the industry's most cost effective plug and play benching systems. It's unique, telescoping frame allows Forte work surfaces to adjust in 1" increments. Height adjustability is complimented with performance; worktops accommodate uniform loads of 1000lb each and upper shelves accommodate 200lbs each.

Single-sided and double-sided configurations are both available—allowing you to select the ideal level of adaptability and value for your specific environment.

Forte's true plug and play service connections ensure that bench and service reconfigurations can always be quickly completed — at minimum cost and down-time.

- 1500 lb. weight load capacity
- Height adjustable worktops: 29"-36" in 1" increments
- Single or double sided configurations
- Integrated plug and play services

As with all A.T. Villa solutions, the Forte benching system stands 100% independently from building structure. UL listed raceways and CSA lab gas services, combined with overhead Ceiling Interface Panels [CIPs] create the ultimate adaptable solution for any lab.



Forte is simplicity and strength at it's best.

MADE IN THE U.S.A.



forte

adaptable laboratory system



Ceiling Interface Panels [CIPs] are configured to the exact requirements of your lab; with the ability to accommodate future expansions.



Caster-based under bench storage provides maximum adaptability and is available in several materials and configurations.

Forte is also available in perimeter, free-standing and sink configurations. A myriad of supporting accessories ensure that your entire laboratory is not only functional but also a true aesthetic showcase.



Widths	36", 48", 60", 72", 84" & 96"
Depths	30" & 60"
Heights	36" & 84"
Worksurfaces	Adjustable 29" - 36": Phenolic Resin, Epoxy Resin, Stainless Steel & Laminate
Pedestals & Overheads	Wood Box/Wood Faces, Steel Box/Steel Faces, Steel Box/Wood Faces & High Density Laminate
Services	UL Listed Raceways, UL Listed Power Strips, CSA Certified Inert Gases & CSA Certified Burning Gases

A.T. VILLA
MADE IN THE U.S.A.


www.atvilla.net
sales@atvilla.net
p 800.554.9259
f 978.582.7488

Contact us today to learn how world-class products, innovative global - design experience and a passion for excellence result in truly superior research environments.