



# SPECIAL INTRODUCTORY PRICING

## SAVE UP TO 35% ON A BUNDLED PACKAGE



A Better Way to Separate

A general purpose bench top centrifuge built to NuAire's high quality standards. NuWind centrifuges offer a great capacity to footprint ratio. With speeds up to **18,000 rpm (30,980 xg)**, a max capacity of 4 x 625 ml, and the ability to store up to 100 programs there isn't an application NuWind can't meet.



## CELL CULTURE PACKAGE


### CHOOSE A CENTRIFUGE:

NuWind NU-C200V - Ventilated



Nothing Over 37° C

NuWind NU-C200R - Refrigerated



Any Accessory at 4° C

### CHOOSE 1 LITER OR 2 LITER APPLICATION PACKAGE:

#### 1 Liter Capacity - C200G214

- Swing out rotor NU-RX250
  - Max speed Refrigerated 4800 rpm/4650 xg
  - Max speed Ventilated 4500 rpm/4650 xg
- 4 x Buckets NU-B250
- 4 x 15ml inserts NU-1T15C (3 x 15ml, 28 Capacity)\*
- 4 x 50ml inserts NU-1T50C (4 x 50ml, 16 Capacity)\*
- Rotor Stand NU-RS2

\* Inserts not compatible with biocontainment lids

#### 2 Liter Capacity - C200G216

- Swing out rotor NU-RX500
  - Max speed Refrigerated 5100 rpm/5350 xg
  - Max speed Ventilated 4500 rpm/4188 xg
- 4 x Buckets NU-B500
- 4 x 15ml inserts NU-2T15C (11 x 15ml, 44 Capacity)\*
- 4 x 50ml inserts NU-2T50C (5 x 50ml, 20 Capacity)\*
- Rotor Stand NU-RS2

\* Inserts not compatible with biocontainment lids



## 30 DAY TRIAL - BEFORE YOU BUY

Visit [www2.nuaire.com/13662](http://www2.nuaire.com/13662) to learn more and request a quote.

Order an additional rotor for a Multifunction centrifuge.

### CLICK\*SPIN

#### ClickSpin Technology:

switch applications within seconds, no tools required.



Your Local NuAire Partner



Call us at (844) 742-3818 or  
Click Here for more options.

