

FIRST YOU CATCH THE WAVE,
THEN YOU GET TO RIDE IT.

INTRODUCING


TWINCORE XLS



-96°C

THE SAFEST ULT FOR HIGH VOLUME AND LONG TERM SAMPLE STORAGE WAS DESIGNED WITH ALL THE FEEDBACK FROM OUR CUSTOMERS IN MIND. IT FEATURES EXTRA-LARGE STORAGE CAPACITY ON A SLIM LINE CONSTRUCTION THAT MAKES IT EASIER TO PLACE IT IN THE LAB.

PULL-DOWN TIME FROM
+20°C to -90°C



CO2
BACKUP SYSTEM
IS NOT ESSENTIAL

7Hr 00 Min to -90
5Hr 24 Min to -82
4Hr 24 Min to -70
Uniformity +/-3°C

UP TO
708
BOX CAPACITY

NEW SLIM LINE DESIGN
THERMAL INSULATED INNER DOORS
TWINCORE TECHNOLOGY
MAXIMUM SAMPLE STORAGE
WIZBOX AND TOUCH SCREEN OR LCD TOUCH CONTROL PANEL
REDESIGNED HANDLE WITH AUTOMATIC VACUUM RELEASE
FINGER PRINT ID LOCK OPTIONAL