

Nuaire PharmaGard ES PR 797 Compounding Pharmacy Isolator

The PharmaGard ES PR797 Compounding Pharmacy Isolator creates a recirculating positive pressure ISO Class 5 clean air controlled environment for aseptic compounding of non-hazardous drugs guaranteed to meet your USP 797 compliance requirements.

Rating: Not Rated Yet
[Ask a question about this product](#)
Manufacturer [Nuaire](#)

Description

-  Ergonomic Design Contributes to a friendly working environment
-  Energy Efficiency Consumes low energy with quality functions
-  Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

- [Specifications](#)

Specifications

Model	NU-PR797-400	NU-PR797-600
Nominal Size	4 ft. (1.2 m)	6 ft. (1.8 m)
Style of Cabinet	Bench top/console with optional base	Bench top/console with optional base stand
Exterior Dimensions (W x D x H)	50 x 32 1/2 x 51 1/2 in. [mm]	74 x 32 1/2 x 51 1/2 in. [mm]
Interior Dimensions (W x D x H)	35 1/4 x 23 1/4 x 27 in. [mm]	59 1/4 x 23 1/4 x 27 in. [mm]
Interchange Dimensions (W x D x H)	14 x 24 x 27 in. [mm]	14 x 24 x 27 in. [mm]

Weight 425 lbs. / 193 kg. 575 lbs. / 261 kg.
Downflow 35 fpm (0.18 m/s) 35 fpm (0.18 m/s)
Velocity
Heat Rejected, 539 880
BTU, Per Hour
Light Intensity 40 to 80 (439 to 40 to 80 (439 to 860)
foot-candles 860)
(LUX)
*Add 3" (76
mm) for base
stand

Warranty 3 Year Parts and
Labor (US & CA)
4 Year Parts (Intl)

Accolades CETA CAG-002-2006, ISO Class 5, UL,
UL-C, CE, Energy Saver, Personnel
Product Environmental Protection, USP
797 Compliant, RoHS
Please note specifications are subject to change. Please
contact NuAire to verify before purchase.

Reviews

There are yet no reviews for this product.