

Nuair PharmaGard ES NR 797 Compounding Containment Pharmacy Isolator

The PharmaGard ES NR797 Compounding Containment Pharmacy Isolator creates negative pressure recirculating ISO Class 5 clean air environment for aseptic compounding of hazardous drugs and minute quantities of chemotherapy drugs to meet USP 797 requirements.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Nuair](#)

Description



Ergonomic Design Contributes to a friendly working environment



Energy Efficiency Consumes low energy with quality function



Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

- [Specifications](#)

Specifications

Model	NU-NR797-400	NU-NR797-600
Nominal Size	4 ft. (1.2 m)	6 ft. (1.8 m)
Style of Cabinet	Bench top/console with optional base stand	Bench top/console with optional base stand
Exterior Dimensions (W x D x H)	50 x 32 ½ x 51 ½ (1270 x 826 x 1308)	74 x 32 ½ x 51 ½ (1880 x 826 x 1308)
Interior Dimensions (W x D x H)	35 ¼ x 23 ¼ x 27 (895 x 591 x 695)	59 ¼ x 23 ¼ x 27 (1505 x 591 x 695)
Interchange	14 x 24 x 27	14 x 24 x 27

Dimensions (W (359 x 610 x 695) (359 x 610 x 695)
x D x H)
in. [mm]
Weight 440 lbs. / 200 kg. 590 lbs. / 268 kg.
Downflow 35 to 45 fpm (.18 35 to 45 fpm (.18 to .23
Velocity to .23 m/s) m/s)
Heat Rejected, 795 1136
BTU, Per Hour
Light Intensity 40 to 80 (439 to 40 to 80 (439 to 860)
foot-candles 860)
(LUX)
*Add 3" (76
mm) for base
stand

Warranty 3 Year Parts and
Labor (US & CA)
4 Year Parts (Intl)

Accolades CETA CAG-002-2006, ISO Class 5, UL,
UL-C, CE, Energy Saver, Personnel
Product Environmental Protection, USP
797 Compliant, RoHS

Please note specifications are subject to change. Please
contact NuAire to verify before purchase.

Reviews

There are yet no reviews for this product.