

Nuaire PharmaGard ES NTE 797 Total Exhaust Compounding Containment Pharmacy Isolator

The PharmaGard ES NTE797 Total Exhaust Compounding Containment Pharmacy Isolator creates negative pressure ISO Class 5 clean air environment for aseptic compounding of hazardous drugs and chemotherapy drugs to meet USP 797 requirements.

Rating: Not Rated Yet
[Ask a question about this product](#)

Manufacturer [Nuaire](#)

Description

-  Ergonomic Design Contributes to a friendly working environment
-  Energy Efficiency Consumes low energy with quality functions
-  Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

- [Specifications](#)

Specifications

Model	NU-NR797-400	NU-NR797-600
Nominal Size	4 ft. (1.2 m)	6 ft. (1.8 m)
Style of Cabinet	Bench top/console with optional base stand	Bench top/console with optional base stand
Exterior Dimensions (W x D x H)	50 x 32 ½ x 51 ½ (1270 x 826 x 1308) in. [mm]	74 x 32 ½ x 51 ½ (1880 x 826 x 1308)
Interior Dimensions (W x D x H)	35 ¼ x 23 ¼ x 27 (895 x 591 x 695) in. [mm]	59 ¼ x 23 ¼ x 27 (1505 x 591 x 695)
Interchange	14 x 24 x 27	14 x 24 x 27

Dimensions (W x D x H)	(359 x 610 x 695)	(359 x 610 x 695)
in. [mm]		
Weight	465 lbs. / 211 kg.	615 lbs. / 279 kg.
Downflow	35 to 45 fpm (.18 to .23 m/s)	35 to 45 fpm (.18 to .23 m/s)
Velocity	.23 m/s)	m/s)
Exhaust	425/722	550/935
CFM/CMH	2.0/w.g./51 mm	2.0/w.g./51 mm w.g.
Standard:	w.g.	
Plan Duct Static Pressure		
Eng/Metric:		
(water gauge)		
Heat Rejected, BTU, Per Hour	539	880
Light Intensity foot-candles (LUX)	40 to 80 (439 to 860)	40 to 80 (439 to 860)
*Add 3" (76 mm) for base stand		
Warranty	3 Year Parts and Labor (US & CA) 4 Year Parts (Intl)	
Accolades	CETA CAG-002-2006, ISO Class 5, UL, UL-C, CE, Energy Saver, Personnel Product Environmental Protection, USP 797 Compliant, RoHS	
Please note specifications are subject to change. Please contact NuAire to verify before purchase.		

Reviews

There are yet no reviews for this product.