

BMT Stericell

111I (3.92 cu ft.) model on rolling cart

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [BMT](#)

Description

- [Description](#)

Description

Discovery Scientific Solutions provides top-of-the-line Stericell lab ovens engineered for hot air sterilization of materials for specific time and temperature ranges.

Equipped with advanced features, our Stericell laboratory ovens operates efficiently and quietly through its patented fine system of forced air convection system.

Patented Forced Air Convection system with temperature range up to 250°C with 300°C option

Hot air Sterilization
Depyrogenation
Drying of components
Drying of non-flammable chemicals
Drying of wet or moist material
160-170-180°C Sterilization cycles

Applications

Research & Laboratories
Pharmaceutical Industry
Animal Science
Industrial Labs
Hospitals

Reviews

There are yet no reviews for this product.