

BMT Venticell IL

1500l (53 cu ft.) model with stainless steel exterior

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [BMT](#)

Description

- [Description](#)

Description

Discovery Scientific Solutions provides BMT Venticell IL lab ovens integrally designed for thorough sterilization from 160°C to 180°C and depyrogenation from 200°C to 300°C temperature ranges of thermally resistant and inflammable materials.

Designed in compliance with the certified quality assurance system EN ISO 9001:2008, Venticell IL is proven to deliver exceptional performance in research and laboratories, pharmaceutical, Bsl 3 / Bsl 4 and industrial laboratories.

Sterilization 160°C to 180°C / Depyrogenation 200°C to 300°C

ISO Class 5 model
ISO Class 7 model
Complies with ISO 14644-1
Hot air Sterilization
Depyrogenation

Applications

Research & Laboratories
Pharmaceutical
Bsl 3 / Bsl 4
Industrial Labs

Reviews

There are yet no reviews for this product.