

Fluorometer Module for DS-11 /DS-11+

Instantly add the power of fluorescence detection to any DeNovix DS-11 or DS-11+ Spectrophotometer.

Rating: Not Rated Yet
[Ask a question about this product](#)

Manufacturer [DeNovix](#)

Description

- [Description](#)
- [Specifications](#)

Description

Add Fluorescence Detection to Any DS-11 or DS-11+

Instantly add the power of fluorescence detection to any DeNovix DS-11 or DS-11+ Spectrophotometer. Simply plug in this USB powered accessory and DeNovix fluorometric EasyApps® are immediately enabled. The easy-to-use apps are preconfigured for use with a wide choice of commercially available assays including those designed for Qubit® and Quantus™ platforms.

- Four Channel Fluorescence
- 0.5 mL thin-walled PCR tubes
- USB Powered

Specifications

FX Module Specifications

Light Sources	Blue LED (~470 nm), Green LED (~525 nm), Red LED (~635 nm)
Excitation Filters	Blue: 442-497 nm, Green: 490-558 nm, Red: 613-662 nm, UV: 361-389 nm
Emission Filters	Channel 1: 514-567 nm, Channel 2: 565-650 nm

Detectors
Tube Type

Channel 3: 664-740 nm, Channel 4: 435-485 nm
Photodiode, detection range 300-1000 nm
0.5 mL Real Time thin-wall PCR tube (polypropylene)

- Four Channel Fluorescence
- EasyApps® with Assay Choices
- Android™ with HD Touchscreen
- Wi-Fi, Ethernet, USB

Reviews

There are yet no reviews for this product.