

dsDNA Fluorescence Quantification Assays

The DeNovix dsDNA Fluorescence Quantification Assays include three highly-specific, highly-sensitive kits for dsDNA quantification. Available in 50, 250 and 1000 assays/kit sizes.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [DeNovix](#)

Description

DeNovix dsDNA Assays

The **DeNovix dsDNA Fluorescence Quantification Assays** include three highly-specific, highly-sensitive kits for dsDNA quantification. Selectively and accurately measure dsDNA in samples at single cell levels even if samples include contaminants such as RNA, ssDNA or protein. Assays are available in 1000 assays/kit size, 250 assays/kit size and 50 assays/kit evaluation size.

- Greatest Dynamic Range 0.5 pg/ μ L – 4000 ng/ μ L
- 20X Greater Sensitivity
- Included in EasyApps® Software
- Two Point Curves, Mix-and-Measure

Reviews

There are yet no reviews for this product.