

BlotBot® 120 Freedom Rocker™

BlotBot 120 Rocker

The BlotBot® 120 is a new Freedom Rocker™ automated blot processor. Using your same protocols, it automates the handling of blots and gels. With microprocessor control, it provides consistent processing conditions. The BlotBot will hybridize/stain and wash your blots/gels during meetings, overnight, or even over a weekend.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Next Advance](#)

Description

The BlotBot® 120 is a new Freedom Rocker™ automated blot processor. Using your same protocols, it automates the handling of blots and gels. With microprocessor control, it provides consistent processing conditions. The BlotBot will hybridize/stain and wash your blots/gels during meetings, overnight, or even over a weekend, so your blots will be ready sooner for downstream analysis. This model provides the same convenience and reliability of previous Freedom Rocker models with the following enhancements:

- Reagent dispensing (e.g. antibody) minimum volume of 3.5 mL
- Holds 1 midi blot, 2 mini blots or strips
- Individual reagent recovery
- Intuitive user interface with color LCD display. You can view and program an entire protocol on one screen.

The Freedom Rocker™ BlotBot® 120 comes with everything you need to get started, including reagent, wash buffer and waste containers. 2 year warranty.

Reviews

There are yet no reviews for this product.