

Hettich Mikro 200 & 200R Microliter Centrifuges

MIKRO 200 & MIKRO 200R tabletop Microliter Centrifuges are ideal for research centrifugation tasks such as PCR, DNA & RNA. The 200R provides fast and highly reliable refrigeration. The Fast Cool Function refrigerates quickly: down to + 4 °C in 10 to 15 minutes, making the MIKRO 200R microcentrifuge ideal for the fast pace laboratory.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Hettich](#)

Description

- [Description](#)
- [Features](#)
- [Specifications](#)
- [Complete Centrifuge Packages](#)

Description

Hettich MIKRO 200 / 200R Microliter Centrifuges

- **FIVE YEAR WARRANTY**
 - All Hettich centrifuges are protected by our 5-year warranty and full technical support.
- **MIKRO 200 APPLICATION PACKAGES**
 - Everything you need for your application in one complete centrifuge bundle. Hettich application packages include the centrifuge, rotor, buckets, adapters and lids.
- **QUIET OPERATION**
 - Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.
- **FUNCTIONAL DESIGN**
 - Easy to lift lid with front-facing control panel for effortless interaction and programming.
- **FULLY CERTIFIED**
 - Engineered and manufactured in accordance with all relevant safety and production standards.
- **ENHANCED SAFETY**
 - Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security
- **MADE IN GERMANY**
 - German engineered and manufactured to be quiet, reliable and safe since 1904.

Everything you need for your application in one complete centrifuge bundle. Hettich application packages include the MIKRO 200 or 200R centrifuge, rotor, buckets, adapters and lids.

- MIKRO 200/200R Microliter Tube Package - 24
- MIKRO 200/200R Microliter Tube Package - 30
- MIKRO 200/200R Microliter Tube Package - Strip
- MIKRO 200/200R CRYO Tube Package
- MIKRO 200/200R Spin Column Package - BC
- MIKRO 200/200R Custom (User Defined) Packages

PACKAGE CONTENTS: Centrifuge, Fixed-angle Rotor, (24) Adapters, Bio-containment Lid

Features

MIKRO 200 / 200R FEATURES

- Robust & steady, temperature control system with pre-cool function from -10 to +40 °C (Optional)
- For DNA, RNA & temperature sensitive samples
- Bio-containment, autoclavable, phenol resistant lids
- Fast cool function: down to +4 °C in 10-15 minutes
- Quiet Operation
- Easy-Lift™ / Safety-Lock™ Lid
- Five (5) year warranty

Specifications

Specifications for MIKRO 200 & 200R Microliter Centrifuges

MIKRO 200 / MIKRO 200R SPECIFICATIONS

Centrifuge	MIKRO 200	MIKRO 200R
Max Capacity	30 x 1.5 / 2.0 mL	30 x 1.5 / 2.0 mL
Max RPM / RCF	15,000 / 21,382	15,000 / 21,382
Temperature Control	n/a	-10 °C to +40 °C
Dimensions (HxWxD)	10.23" x 11" x 13.5"	10.23" x 11" x 21.53"
Weight	25.35 lbs	61.72 lbs
Cat. No. 115V	2400-01	2405-01

Additional Hettich Options: MIKRO 185/220/220R

Centrifuge	MIKRO 185	MIKRO 220 / 220R
Max Capacity	24 x 1.5 / 2.0 mL	48 x 0.2 - 2.0 mL
Max RPM / RCF	13,300 / 17,008	18,000 / 31,514
Temperature Control	n/a	n/a / -20 °C to +40 °C
Dimensions (HxWxD)	9" x 10.3" x 13.9"	12.3" x 13" x 16.5" / 12.3" x 13" x 25.5"
Weight	22 lbs	46.29 lbs / 95.59 lbs
Cat. No. 115V	1203-01	2200-01 / 2205-01

Complete Centrifuge Packages

MIKRO 200 & 200R Centrifuge Packages for Common Applications

Everything you need for your application bundled into one complete centrifuge package, ready to go. Hettich application packages include MIKRO 200 or 200R centrifuge, rotor, buckets, adapters and bio-containment lids (when applicable).

MIKRO 200/200R Microliter Tube Package - 24

TUBE QTY	VOLUME (mL)	SIZE (mm)	RPM	g-Force
24	1.5 mL - 2.0 mL	11 x 38 mm	15,000 RPM	21,382 RCF

PACKAGE CONTENTS:

(1) Centrifuge, (1) Fixed-angle Rotor, (24) Adapters 1.5 mL
Add a Bio-containment Lid (optional).
0.2mL – 0.8 mL adapter combinations are available.

MIKRO 200/200R Microliter Tube Package - 30

TUBE QTY	VOLUME (mL)	SIZE (mm)	RPM	g-Force
30	1.5mL - 2.0 mL	11 x 38 mm	15,000 RPM	20,627 RCF

PACKAGE CONTENTS:

(1) Centrifuge, (1) Fixed-angle Rotor, (30) Adapters 1.5 mL
Add a Bio-containment Lid (optional).
0.2mL – 0.8 mL adapter combinations are available.

MIKRO 200/200R Microliter Tube Package - Strip

TUBE QTY	VOLUME (mL)	SIZE (mm)	RPM	g-Force
4	PCR Strip (8)	--	15,000 RPM	14,338 RCF

PACKAGE CONTENTS:

(1) Centrifuge, (1) Fixed-angle Strip Rotor

MIKRO 200/200R CRYO Tube Package

TUBE QTY	VOLUME (mL)	SIZE (mm)	RPM	g-Force
20	1.8 mL	--	15,000 RPM	21,130 RCF

PACKAGE CONTENTS:

(1) Centrifuge, (1) Fixed-angle Rotor
Add a Bio-containment Lid (optional).

MIKRO 200/ 200R Spin Column Package - BC (Biocontainment)

TUBE QTY	VOLUME (mL)	SIZE (mm)	RPM	g-Force
24	1.5 - 2.0 mL	11 x 38 mm	15,000 RPM	21,382 RCF

PACKAGE CONTENTS:

Centrifuge, Fixed-angle Rotor, (24) Adapters, Biocontainment Lid

Need a Different Centrifuge Package? - Let Us Know!

Create Your Own Custom Package

Additional packages and configurations are available. We offer a complete centrifuge product range from small microliter centrifuges, compact tabletop centrifuges, medium to large benchtop centrifuges, high capacity benchtop and floor standing centrifuges for general purpose and specialty applications. Submit your application specific requirements using the "Request Quote" form above for custom package details to meet your needs.

Reviews

There are yet no reviews for this product.