

Cardboard Boxes for LN2 & ULT Freezer Storage Racks

Water repellent cardboard boxes for storage of samples or for tubes are supplied with or without cell dividers that have drain holes in each corner.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Discovery Scientific Solutions](#)

Description

Our white or yellow water repellent cardboard boxes for tube storage are supplied with or without cell dividers. Dividers are available with 100, 81, 64, 49, 36, 25 or 16 slots. These standard boxes have four drain holes in each corner.

The white water repellent cardboard boxes for 15ml or 50ml centrifuge tubes are supplied with or without 16 or 36-slot Cell dividers. These boxes also have four drain holes in each corner.

The white water repellent mini-cardboard boxes for storage of samples or tubes are supplied with or without 25-Cell dividers. Mini-boxes also have drain holes.

Part Number	Description	Dimension L/D/H
SB2CC	Standard White 2" Cardboard Box	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm
SB2CC-XX*	Standard White 2" Box w/cell divider	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm
SB2C	Standard Yellow 2" Cardboard Box	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm
SB2C-XX*	Standard Yellow 2" Cardboard Box w/cell divider	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm
SB3C	Standard White 3" Cardboard Box	5 1/4 x 5 1/4 x 2 15/16 inch 133 x 133 x 75mm
SB3C-XX*	Standard White 3" Cardboard Box w/cell divider	5 1/4 x 5 1/4 x 2 15/16 inch 133 x 133 x 75mm
SB4C	Standard White 4 3/4" Cardboard Box with Drain Hole	5 3/4 x 5 3/4 x 4 3/4 inch 146 x 146 x 121mm
SB4C-XXX**	Standard White 4 3/4" Cardboard Box w/Cell Divider & Drain Hole	5 3/4 x 5 3/4 x 4 3/4 inch 146 x 146 x 121mm
MB2C	White 2" Mini Cardboard Box with Drain Slots	3 x 3 x 2 inch 75 x 75 x 51mm
MB2C-25	White 2" Cardboard Mini Box w/Cell Divider & Drain Slots	3 x 3 x 2 inch 75 x 75 x 51mm
MB3C	White 3" Mini Cardboard Box with Drain Slots	3 x 3 x 3 inch

Part Number	Description	Dimension L/D/H
MB3C-25		75 x 75 x 76mm
	White 3" Mini Cardboard Box w/Cell Divider &	3 x 3 x 3 inch
	Drain Slots	75 x 75 x 76mm

Reviews

There are yet no reviews for this product.