

Benchmark MC-24 Touch Microcentrifuge

Compact design with full color touch screen. High speed (up to 13,500 rpm/16,800 xg) with 1 button program recall. Combi-Rotor accepts tubes and PCR strips, making this the perfect partner for any molecular biology laboratory. 2 year warranty.

Rating: Not Rated Yet
[Ask a question about this product](#)
Manufacturer [Benchmark Scientific](#)

Description

MC-24 Touch Microcentrifuge

Benchmark's MC-24 Touch is the first microcentrifuge in its class to feature touch screen control. In addition to the next generation control panel, the centrifuge also offers high speed (up to 13,500 rpm/16,800 xg) and a unique COMBI-Rotor™ that accepts both standard microtubes (1.5/2.0ml) and PCR tubes/strips (0.2ml), making this the perfect partner for any molecular biology laboratory.

The large, color touch screen shows actual speed (in rpm or rcf) along with set and remaining time. Up to 9 speed/time profiles can be stored in memory and the home screen displays a "mySpin" icon for instant recalls a designated favorite profile. Supplied with Benchmark's unique COMBI-Rotor™, the MC-24 Touch will accept 24 x 1.5/2.0ml tubes and 2 strips of 0.2ml tubes (16 x 0.2ml) without requiring separate adapters. Samples are kept cool, even during extended runs at the maximum speed setting by a proprietary air flow system that moves air evenly and quietly throughout the chamber. The MC-24 Touch is cold room safe and carries a 2 year warranty.

- Full Color Touch Screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact design, 9 x 12 x 7.75 in.
- Unique air-flow system, keeps samples cool

Specifications

Speed	0,000 to 13,500 rpm (100 to 16,800 xg)
Capacity	24 x 1.5/2.0ml
Increment	2 x PCR Strips (16 x 0.2ml)
Timer	100 rpm/ 100 rcf
Momentary	1 to 99 min./cont.
Dimensions	Yrs
Weight	8 x 12 x 7.75 inches
Electrical	14.5lb / 6.6kg
Warranty	115V, 60 Hz or 230V, 50 Hz
	2 Year

Reviews

There are yet no reviews for this product.