

## Helmer iLR104-ADA Undercounter Laboratory/Pharmacy Refrigerator

The iLR104-ADA is a single door ADA compatible Undercounter i.Series® Laboratory/Pharmacy Refrigerator. It fits under 34" ADA compliant counter heights. It features the i.C<sup>3</sup> full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

### Description

- [Description](#)
- [Features](#)
- [Specification](#)

## Description

The iLR104-ADA is a single door ADA compatible Undercounter i.Series® Laboratory/Pharmacy Refrigerator. It fits under 34" ADA compliant counter heights. It features the i.C<sup>3</sup> full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

## Features

- Innovative i.C<sup>3</sup> User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Rel.i Plus warranty includes 7 years compressor, 2 years parts, and 1 year labor

## Specification

### Interior Dimensions (w x h x d)

- 19.5" x 17.75" x 22"
- 19.5" x 17.75" x 22"

### Exterior Dimensions (w x h x d)

- 24" x 31.5" x 26.75"
- 610mm x 801mm x 680mm

### Interior Volume (w x h x d)

- 4 cu ft / 113 liters

### Shelves (w x d)

- 2 Adjustable Full-Size
- 18" x 18" (458mm x 458mm)

### Electrical

- 115V 60Hz 5 FLA

### Net Weight

- 191 lb / 87 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

## Reviews

There are yet no reviews for this product.