

## Helmer HLR245 Double Door Laboratory Refrigerator

The HLR245 is a double door Horizon Series Laboratory Refrigerator. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

### Description

The HLR245 is a double door Horizon Series Laboratory Refrigerator. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

### Features

- Temperature controller with alarm/monitor system provides protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i warranty includes 5 years compressor, 2 years parts, and 1 year labor

## Specification

### Interior Dimensions (w x h x d)

- 55" x 58.25" x 24.25"
- 1397mm x 1480mm x 616mm

### Exterior Dimensions (w x h x d)

- 59.25" x 78.75" x 29.5"
- 1505mm x 2001mm x 750mm

### Interior Volume (w x h x d)

- 44.9 cu ft / 1271 liters

### Shelves (w x d)

- 8 Adjustable Full-Size
- 22" x 23" (559mm x 585mm)

### Electrical

- 115V 60Hz 11.5 FLA
- 230V 50Hz 6 FLA
- 230V 60Hz 6 FLA

### Net Weight

- 701 lb / 318 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

## Reviews

There are yet no reviews for this product.