

Helmer iLF120 Laboratory Freezer

The iLF120 is a single door i.Series® Laboratory Freezer. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The iLF120 is a single door i.Series® Laboratory Freezer. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Features

- Innovative i.C³ User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Rel.i warranty includes 5 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 24.75" x 58.25" x 24.25"
- 629mm x 1480mm x 616mm

Single Door Upright: Helmer iLF120 Laboratory Freezer

Exterior Dimensions (w x h x d)

- 29.5" x 80" x 29.5"
- 750mm x 2032mm x 750mm

Interior Volume (w x h x d)

- 20.2 cu ft / 572 liters

Shelves (w x d)

- 4 Adjustable Full-Size
- 22" x 23" (559mm x 585mm)

Electrical

- 208 / 230V 60Hz 4.3 FLA
- 208 / 230V 60Hz 4.3 FLA

Net Weight

- 443 lb / 201 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

Reviews

There are yet no reviews for this product.