

Helmer iLF125 Laboratory Freezer

The iLF125 is a single door i.Series® Laboratory Freezer. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The iLF125 is a single door i.Series® Laboratory Freezer. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Features

- Innovative i.C³ User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Rel.i warranty includes 5 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 24.75" x 58.25" x 30.25"
- 629mm x 1480mm x 769mm

Single Door Upright: Helmer iLF125 Laboratory Freezer

Exterior Dimensions (w x h x d)

- 29.5" x 80" x 35.5"
- 750mm x 2032mm x 902mm

Interior Volume (w x h x d)

- 25.2 cu ft / 714 liters

Shelves (w x d)

- 4 Adjustable Full-Size
- 22" x 29" (559mm x 737mm)

Electrical

- 208 / 230V 60Hz 4.3 FLA
- 230V 50Hz 3.8 FLA

Net Weight

- 481 lb / 219 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

Reviews

There are yet no reviews for this product.