

Helmer iPR256 Double Door Pharmacy Refrigerator

The iPR256 is a double door i.Series® Pharmacy Refrigerator. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The iPR256 is a double door i.Series® Pharmacy Refrigerator. It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Features

- Innovative i.C³ User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i Plus warranty includes 7 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 55" x 58.25" x 30.25"

Double Door Upright: Helmer iPR256 Double Door Pharmacy Refrigerator

- 1397mm x 1480mm x 769mm

Exterior Dimensions (w x h x d)

- 59.25" x 79.5" x 35.5"
- 1505mm x 2020mm x 902mm

Interior Volume (w x h x d)

- 56 cu ft / 1586 liters

Shelves / Baskets (w x d)

- 2 Adjustable Full-Size
- 22" x 29" (559mm x 737mm)
- 12 Roll Out
- 21" x 27.25" (534mm x 693mm)

Electrical

- 115V 60Hz 11.5 FLA
- 230V 50Hz 6 FLA
- 230V 60Hz 6 FLA

Net Weight

- 876 lb / 398 kg

Certified to applicable UL and CSA standards by a NRTL

- CE 0086

Reviews

There are yet no reviews for this product.