

Helmer iPR225 Pass-Thru Pharmacy Refrigerator

The iPR225 is a two door i.Series® Pass-Thru Pharmacy Refrigerator (one door on each side). It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The iPR225 is a two door i.Series® Pass-Thru Pharmacy Refrigerator (one door on each side). It features the i.C³ full-color touchscreen monitor with downloading, maximum temperature uniformity, and superior cabinet construction.

Features

- Innovative i.C³ User Interface and monitor offers advanced monitoring for protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i Plus warranty includes 7 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 24.75" x 58.25" x 32"

Double Door Upright: Helmer iPR225 Pass-Thru Pharmacy Refrigerator

- 629mm x 1480mm x 813mm

Exterior Dimensions (w x h x d)

- 29.5" x 80" x 36"
- 750mm x 2032mm x 915mm

Interior Volume (w x h x d)

- 26.5 cu ft / 750 liters

Shelves / Baskets (w x d)

- 2 Adjustable Full-Size
- 22" x 29" (559mm x 737mm)
- 3 Solid Stainless Steel Two-way
- 22" x 27" (559mm x 686mm)

Electrical

- 115V 60Hz 9.25 FLA

Net Weight

- 534 lb / 243 kg

Certified to applicable UL and CSA standards by a NRTL

Reviews

There are yet no reviews for this product.