

Helmer HPR245 Double Door Pharmacy Refrigerator

The HPR245 is a double door Horizon Series Pharmacy Refrigerator. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Helmer](#)

Description

- [Description](#)
- [Features](#)
- [Specification](#)

Description

The HPR245 is a double door Horizon Series Pharmacy Refrigerator. It offers a temperature controller with alarm/monitor, maximum temperature uniformity, and superior cabinet construction.

Features

- Temperature controller with alarm/monitor system provides protection of stored products
- Powerful refrigeration system provides maximum temperature uniformity and quick recovery
- Superior cabinet built with high-quality construction meets or exceeds performance standards required by the medical and scientific industries
- Energy-efficient LED lighting spotlights storage providing maximum illumination
- Rel.i warranty includes 5 years compressor, 2 years parts, and 1 year labor

Specification

Interior Dimensions (w x h x d)

- 55" x 58.25" x 24.25"

Double Door Upright: Helmer HPR245 Double Door Pharmacy Refrigerator

- 1397mm x 1480mm x 616mm

Exterior Dimensions (w x h x d)

- 59.25" x 78.75" x 29.5"
- 1505mm x 2001mm x 750mm

Interior Volume (w x h x d)

- 44.9 cu ft / 1271 liters

Shelves / Baskets (w x d)

- 2 Adjustable Full-Size
- 22" x 23" (559mm x 585mm)
- 12 Roll Out
- 21" x 21.25" (534mm x 540mm)

Electrical

- 115V 60Hz 11.5 FLA
- 230V 50Hz 6 FLA
- 230V 60Hz 6 FLA

Net Weight

- 823 lb / 374 kg

Certified to applicable UL and CSA standards by a NRTL

- CE0086

Reviews

There are yet no reviews for this product.