

## Steelco AC 1300

The device is an efficient system able to clean and disinfect all types of feeding bottles, animal and fish cages.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Steelco](#)

### Description

- [Description](#)
- [Specification](#)

## Description

Steelco CBW equipment cleans and disinfects all types of feeding bottles, animal and fish cages, assuring outstanding washing and sanitization results.

Final rinse temperature is guaranteed to reach 85°C/185°F through an independent rinsing circuit.

Steelco AC 1300 encloses dedicated spray jets for the bottles and self-cleaning filter for the recirculation of the washing water.

- Separated washing and rinsing circuits.
- Easy access with sliding electric control panel.
- Detergent cabinet with capacity up to two 10 lt. / 2.64 US gal containers.

## Specification

Washers dimensions: width 1760 mm / 69.29"; depth 830 mm / 32.68"; height 1900 mm / 74.80"

Washing chamber capacity: width 1550 mm / 61.02"; depth 730 mm / 28.74"; height 1130 mm / 44.49"

- Drop down/vertical sliding loading door.
- Chamber and doors entirely made in stainless steel AISI 316L, external body made in stainless steel AISI 304.
- Rotating washing and rising arms positioned on the top, middle and bottom levels of the washing chamber.
- Multi stage filtering system.
- Washing pump pressure switch.

- Tank and sump temperature check.
- Equipped with 1 dosing pump of chemical product, controlled by a micro-processor which allows to modify the required product quantity for each type of cycle.
- Check of the chemical product quantity through a minimum level alarm.
- RS 232 port for printer connection.

### Reviews

There are yet no reviews for this product.