

Nuaire AllerGard NU-619 Innovive Edition ATS for INNOCAGE

The AllerGard NU-619 Innovive Edition ATS for INNOCAGE® mouse allows technicians to prepare the work zone under ISO Class 4 conditions. A stack of 25 INNOCAGE® MOUSE bottoms at 13" (330 mm) tall fits neatly into deep wells provided by the NU-619.

Rating: Not Rated Yet

[Ask a question about this product](#)

Manufacturer [Nuaire](#)

Description



Ergonomic Design Contributes to a friendly working environment



Energy Efficiency Consumes low energy with quality functions



Dynamic Air Block Assures strong inflow air velocity and maximum personnel protection

Reviews

There are yet no reviews for this product.